

# EZe-Band® LF

Latex-free self-closure bandage that secures with front and end closures.



# EZe-Band® LF

## Latex-Free Self-Closure Elastic Bandage

EZe-Band® LF is a latex-free, orthopedic grade knitted bandage that wraps and secures easily. It provides patient comfort and is durable enough to be used again and again. EZe-Band LF's front closure anchors the bandage on the first wrap quickly and easily. Its end closure conveniently secures the bandage at any given point without pins, clips or tape. It offers latex-sensitive patients and health professionals an alternative to latex rubber elastic bandages to help prevent possible future latex sensitization.



*EZe-Band® LF's end closures safely secure the bandage without pins, clips or tape.*

### Description

- ▶ Latex-free
- ▶ Double closure delivers maximum control and consistent compression
- ▶ Requires no pins, clips or tape
- ▶ Unique knitted construction provides superior absorbency and breathability
- ▶ Sterility assured by Tyvek® pouch

### Applications

- For controlled and consistent compression and support
- For use in Hospitals, Physician Offices, Industrial Health, Sports Medicine, Long Term Care, Home Health, School Health and Burn Units

#### EZe-Band® LF

Order No.	Width	Stretched Length	Bandages Per Pack	Bandages Per Case
59120000	2"	5 yds.	10	60
59130000	3"	5 yds.	10	60
59140000	4"	5.5 yds.	10	60
59160000	6"	5.5 yds.	10	60

#### Double Length

Order No.	Width	Stretched Length	Bandages Per Pack	Bandages Per Case
59180000	4"	11 yds.	12	36
59190000	6"	11 yds.	12	24

#### Sterile

Order No.	Width	Stretched Length	Bandages Per Pack	Bandages Per Case
59720000	2"	5 yds.	1	36
59730000	3"	5 yds.	1	36
59740000	4"	5.5 yds.	1	36
59760000	6"	5.5 yds.	1	36

#### Sterile Double Length

Order No.	Width	Stretched Length	Bandages Per Pack	Bandages Per Case
59790000	6"	11 yds.	1	24